

Medtronic

Sprinter™ Legend™ 1.25 mm
Balloon Dilatation Catheter

Designed to cross

Sprinter Legend 1.25 mm
balloon features **Zerofold
technology** that:

- Enables a 0.5 mm crossing profile¹
- Has no shoulders or wrapped material for excellent recross

Crossing profile =
0.5 mm/0.020 in¹

Lesion entry profile =
0.41 mm/0.016 in²

No balloon
shoulders



Sprinter Legend
1.25 mm balloon
with Zerofold
technology



Sprinter Legend
3.0 mm balloon with
Miniwrap folding

Sprinter Legend RX balloon ordering information

Balloon diameter (mm)	Balloon length (mm)					Crossing profile O.D. mm (in)
	6	10	12	15	20	
1.25	SPL12506X	SPL12510X	SPL12512X	SPL12515X	SPL12520X	0.5 (0.020)

Sprinter Legend OTW balloon ordering information

Balloon diameter (mm)	Balloon length (mm)					Crossing profile O.D. mm (in)
	6	10	12	15	20	
1.25	SPL12506WL	SPL12510WL	SPL12512WL	SPL12515WL	SPL12520WL	0.5 (0.020)

Sprinter Legend RX/OTW balloon compliance information

Pressure (atm)	1.25 mm balloon diameter
4	0.60
5	0.60
6	0.63
7	0.67
8	0.75
9	0.86
10	0.99
11	1.11
12	1.20
13	1.28
14	1.35
15	1.40
16	—



Nominal pressure and rated burst pressure are the same. Do not exceed.

Nominal pressure: the pressure at which the balloon reaches its labeled diameter

Rated burst pressure: the maximum pressure to which a balloon is designed to be inflated

¹ TR-DA08-020 Competitor Evaluation Report for Sprinter Legend 1.25 OTW balloon dilatation catheter, March 2, 2009.

² 10880534DOC Competitor Evaluation Test Report for the Euphora™ RX balloon dilatation catheter, January 25, 2019.

Important Information:

Prior to use, refer to the product Instructions for Use supplied with these devices for indications, contraindications, side effects, suggestion procedure, warnings, and precautions.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a licensed healthcare practitioner. See package inserts for full product information.

Medtronic

Medtronic
Tel: 707.525.0111

LifeLine Customer Support
Tel: 877.526.7890

Product Services
Tel: 888.283.7868

[medtronic.com](https://www.medtronic.com)